

Device Schedule:

Intended Purpose as per the Instructions for Use:

Coated VICRYL™ Plus Antibacterial Suture is indicated for use in general soft tissue approximation and/or ligation, except for ophthalmic, cardiovascular and neurosurgical tissues

Device Name	Model	Type (Codes as per (EU) 2017/2185)	Risk Classification	Basic UDI-DI
Vicryl Plus Antibacterial Polyglactin 910 Braided Absorbable Suture	Vicryl Plus Antibacterial	MDN 1104	Class III, Implantable	0705031a007925K

Additional Information:

Surgical needle and suture combinations from within the following limits define the device range for the Coated VICRYL™ Plus Antibacterial (Polyglactin 910) Sterile Synthetic Absorbable Suture.

Suture Characteristics	Range
Suture Material (Absorbable/Non-Absorbable)	Absorbable
Suture Gauge Size	1.0 – 5.0 (Metric)
Suture Length	5 cm – 250 cm
Suture Dyed/Undyed	Dyed/Undyed
Suture Color (If dyed)	Violet
Coated/Uncoated	Coated (Copolymer of glycolide and lactide, calcium stearate)
Contains Antimicrobials (Yes/No)	Yes
Triclosan Maximum Levels (µg/m)	≤ 275 µg/m
Multifilament/Monofilament	Multifilament
Accessories to suture type	N/A

Suture Characteristics	Range
Needled/Non-Needled	Needled/Non-Needled
Number of Needles per Suture	Single Armed/Double Armed
Needle Material	420 SS, 4310 SS and ETHALLOY
Needle Coating	Silicone, MULTIPASS, CERBERUS
Needle Shape	Straight/Curve
Needle Color	Silver/Black
Needle Length	3.5 mm – 110 mm
Needle Wire Diameter	0.10 mm – 1.55 mm